

2nd Generation


Natural Biomedical Wound Dressing Solution

Available size: 5 X 5 , 10 X 10

Alginate Wound Dressing

non-woven type

A soft and integral gel is formed when this dressing contacts wound exudate. Excess exudate would be absorbed, and the healing process is significantly improved.



Features

Accelerates wound healing

Promotes cell growth, skin generation, and reduces scarring.

Remarkable absorbency & easy removal

Allows the dressing to remain for a couple of days by absorbing several times its own weight in exudate, and can be easily removed or washed out.

Conforms to wound contours

Soft, comfortable, and easy to be cut for a variety of wound sizes.